

Biopsy Punch



Biopsy punch has a sharp stainless steel cutting edge that is seamless.

Handle has ribs for a secure grip and control.

Each biopsy punch has a size embossed on it for easy identification.

This product is sterile, disposable, ready to use with a tip cover.

Curette Punch



A dermal curette is a spoon-shaped instrument with a sharp inside and a blunt outside.

It is used to scrape and remove skin lesions, growths, and cancers.

Each curette punch has a size embossed on it for easy identification.

This product is stainless steel, sterile, disposable, ready to use with tip cover.

Plunger Biopsy Punch



The benefit of these special punches is that they guarantee a simple way of ejecting lodged tissue using the plunger system.

Handle has ribs for a secure grip and control.

The punch's plunger system is similar to a ball pen.



ALL OUR PRODUCTS ARE STERILE GAMMA UNLESS PACKAGE IS DAMAGED OR OPENED, DISPOSABLE, SINGLE USE ONLY, LATEX-FREE. STORE AT TEMPERATURE BELOW 30° C (80F). THEY ARE PACKED IN A DISPLAY INNER BOX OF 25 UNITS EACH, IN INDIVIDUALLY STERILE POUCHES.

COLOUR CODED	PRODUCT DESCRIPTION	PACK SIZE
Biopsy Punch		
D-BP2M	Stainless Steel Biopsy Punch Individually Wrapped Sterile 2mm with Plastic Handle Single Use	25 pc/box
D-BP3M	Stainless Steel Biopsy Punch Individually Wrapped Sterile 3mm with Plastic Handle Single Use	25 pc/box
D-BP4M	Stainless Steel Biopsy Punch Individually Wrapped Sterile 4mm with Plastic Handle Single Use	25 pc/box
D-BP5M	Stainless Steel Biopsy Punch Individually Wrapped Sterile 5mm with Plastic Handle Single Use	25 pc/box
D-BP6M	Stainless Steel Biopsy Punch Individually Wrapped Sterile 6mm with Plastic Handle Single Use	25 pc/box
Plunger Biopsy Punch		
D-BPP1.5M	Stainless Steel Biopsy Punch with Plunger Individually Wrapped Sterile 1.5mm with Plastic Handle Single Use	25 pc/box
D-BPP2M	Stainless Steel Biopsy Punch with Plunger Individually Wrapped Sterile 2mm with Plastic Handle Single Use	25 pc/box
D-BPP3M	Stainless Steel Biopsy Punch with Plunger Individually Wrapped Sterile 3mm with Plastic Handle Single Use	25 pc/box
D-BPP4M	Stainless Steel Biopsy Punch with Plunger Individually Wrapped Sterile 4mm with Plastic Handle Single Use	25 pc/box
Curette Punch		
D-DC2	Stainless Steel Dermal Curettes 2mm with Plastic Handle Single Use	25 pc/box
D-DC3	Stainless Steel Dermal Curettes 3mm with Plastic Handle Single Use	25 pc/box
D-DC4	Stainless Steel Dermal Curettes 4mm with Plastic Handle Single Use	25 pc/box